

Exhibit 3

BEST AVAILABLE COPY

2/2 Digest FT7 DNA 3, 10, 14 w/ Kpn I / Hind III
use previous prep as control (-)

Expected Size

(3) Sense (4) Antis. (-)

3159	3159
1573	1366
625	782

10	3159	3159
	1573	1124
	650	832
	380	650

14	3159	3159
	1573	945
	650	832
	204	650

Major problem w/ #14

Some how samples got mixed-up

60 band is checked out - Sheet for single colonies

A B : E

Also mini prep from original culture

Digest w/ Kpn I / Hind III 0/14

4/23 Run gel of Digests



Clearly Sample A which was given as mate prep is in the wrong orientation

Start 0/14 of D : E to mini prep before start of 500ml culture

That mini prep FT7 DNA 14 - D : E

Digest w/ Kpn I / Hind III

run on gel - Both are fine
use -D for mate




BEST AVAILABLE COPY

1/25 Plasmid isolation of FT7, cDNA 14 - D up to Binding
 - do 500 nls in 2 250 sets, extract 2 tubes to Band
 7/26 Pull Bands - double band 1 pip - 6 hrs during the Day

7/27 Digest FT7, cDNA 14 w/ Kpn I / Hma III

- ① - Single bandy
- ② Double Bandy
- ③ FT7, cDNA 14 antisense



Check Absorbance / Concentration of


Samples	260	FT7	cDNA 14	260/280
①	.127	.062	.63 ng/ul	2.0
②	.143	.065	.71 ng/ul	2.3

Sequence
 FT7, cDNA 3
 T7

cDNA 14
 T7
 8850
 9007
 8874

Digest Mouse FT4

	260	280	260/280
psr7 genomic	.085	.045	
with Kpn I Coding	.106	.060	
Digest of cDNA 3			
Kpn I cDNA 10			
Hma III cDNA 14			



BEST AVAILABLE COPY

8/4 FACS Analysis of Transfected cells w/ following Vectors

pCDNA 7

FIT7, 1, 2a, 2b, 3

1, 2b, 3

1, 3

cDNA 3

cDNA 10

cDNA 14

2 plates / Vector - Divide

FACS 7mbs

yuko 5mbs

EAT 3mbs

ETassy 5mbs

(75ul)

(100ul)

Antibodies

IgM

H

- blank

1:100

IgM

hxf

green

1:1000

IgM

SLx

green

1:200

IgG

ha

red

1:500

IgG

SLa

blue

1:500

2nd Antibody

IgM 2.5mbs

12.5 / 2.5mbs

IgG 1.5

6.0 / 1.5mbs

Results are

H - all neg

hxf - all neg

SLx, pD (-), 1, 2a, 2b, 3(+), 1, 2b, 3(+), 1, 3 (-), cDNA 3(-), 10 (+), 14 (+)

ha - all neg

SLa - all neg

9/2 Spin 12 punn KG, FT4

6451

6080

2470

6199

6374

6087

6306

6086

6203

6085

5721

5671

9/9 Run sig gel of above samples

Also spin

5728

6084

7213

5731

5737

6082

5662

5727

6201

5725

6200

BEST AVAILABLE COPY

8/10 Spinning gel of 8/9 samples
Sequence

FT4 6079

6202

6307

6373

CDNA10 715

946

CDNA14

T7

8850

8807

8874

Protein assay of FACS Samples, Also CAT assay

pcDNA

FT7 1,2a2b,3

FT7 1,2b,3

FT7 1,3

CDNA3

CDNA10

CDNA14

BSA Blank Protein

0 1.00 116 108

1 0.208 225 214

2 0.369 383 376

4 0.691 673 682

8 1.215 1.230 1.222

16

Sample

pcDNA .292 .294 .293

1,2a2b,3 .337 .330 .33

1,2b,3 .298 .343 .32

1,3 .369 .372 .37

CDNA3 .363 .379 .37

CDNA10 .306 .298 .3

CDNA14 .225 .253 .22

FACS Results

Only stain w/ 5h4

1,2a2b,3 23.6%

1,2b,3 24.6%

CDNA10 14.9%

CDNA14 8.0%

Micro BCA Protein Assay

Reagent mic	MC	MB	MA
Per assay tube (ml)	0.01	0.24	0.25
Cocktail for	Tubes		

Incubate 1 h at 60°C and cool to room temp.

Since the color development has no end point, all tubes must be heated and cooled at the same time

1 mg/ml BSA (l)	Water (l)	Reagent (l)	Abs. 562	
0.0	500.0	500.0	Blank	
1.0	499.0	500.0	0.108	Slope = 0.0734
2.0	498.0	500.0	0.214	Y intercept = -0.0656
4.0	496.0	500.0	0.376	X intercept = -0.8940
8.0	492.0	500.0	0.682	R = 0.9985
16.0	484.0	500.0	1.222	

8/10 Spinning of FT4 9/10 Samples

Sample	l in assay	Water (l)	Reagent (l)	Abs. 562	mg protein/ml
pcDNA1	5.00	495.00	500	0.293	0.62
FT7 1,2a2b,3	5.00	495.00	500	0.333	0.73
FT7 1,2b,3	5.00	495.00	500	0.320	0.69
FT7 1,3	5.00	495.00	500	0.370	0.83
CDNA 3	5.00	495.00	500	0.371	0.83
CDNA 10	5.00	495.00	500	0.302	0.64
CDNA 14	5.00	495.00	500	0.239	0.47

BEST AVAILABLE COPY

9. CAT Assay (FT7) samples

2.5 μ l of cell extract

Control

pcDNA

1,2a2b,3

1,2b,3

1,3

CDNA 3

CDNA 10

CDNA 14

Cocktail 15

3/4 Chlase 300

Tris, 2M pH 8.0 75

But Co A 75

H₂O 300

50 / tube

CAT

Vector		Counts/5ul 8/12/94		Incorporated Counts (.95)	Total Counts		Total Counts Incorporated	
pcDNA1		11,349	11,829	9,189	8,083	238,169	244,643	9,673
FT7 1,2a2b,3		11,161	11,441	27,211	21,919	250,431	250,739	28,643
FT7 1,2b,3		11,772	11,826	37,541	40,684	272,981	277,204	39,517
FT7 1,3		11,215	11,690	23,076	28,706	247,376	262,506	24,291
cdNA 3		11,834	11,206	33,885	39,098	270,565	263,216	35,668
cdNA 10		12,017	11,312	30,066	33,165	270,406	259,405	31,848
cdNA14		11,079	11,570	44,133	40,529	265,713	271,929	46,456
Control			10,354		424		207,504	446
	Protein Conc. (ug/2.5ul)	Total Counts Inc- Bkg		% INC/hr		% INC/hr/ug		Mean CAT Activity
pcDNA1	1.55	9,249	8,063	3.92	3.30	2.53	2.13	2.33
FT7 1,2a2b,3	1.83	28,219	22,649	11.27	9.03	6.17	4.95	5.56
FT7 1,2b,3	1.72	39,093	42,401	14.32	15.30	8.33	8.89	8.61
FT7 1,3	2.07	23,867	29,793	9.65	11.35	4.66	5.48	5.07
cdNA 3	2.07	35,244	40,730	13.03	15.47	6.29	7.48	6.88
cdNA 10	1.60	31,224	34,487	11.55	13.29	7.22	8.31	7.76
cdNA14	1.17	46,032	42,238	17.32	15.53	14.81	13.28	14.04

9. Assemble Data of FT7 FACS, CAT Assay / give to Judy
Work on FT4 Sequencer

9.5 7 deaza sequencing on Troublesome FT4 samples

6451 6200

6378 6079

6306 6307

6203 1897

5721 1898

1899

9.6 Sequencing gel of 8/15 samples (FT7) Formamide gel
Prove Deonias' rather than But w/ GAP prove
To check condition of RFL

9.7 The 7 deaza technique didn't resolve all of the compressions
Try a Terminal transferase technique
Run standard Syntexase rxn, after extension made
Heat tubes (A, C, G, T) for 1.5 mins 100°C
Hold on ice 10 min, Purpore TdT/dNTP cocktail
Add to tubes, 37°C 30 min
Add Spl Stop